TMX[®] META

by BlackBoard[®]

KNOW-HOW TRIFFT KNOW-HOW

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
by BlackBoard[®]

KNOW-HOW MEETS KNOW-HOW


In cooperation with **BlackBoard[®]**, we have succeeded in creating a unique product which helps you to take matters concerning foot into your own hands.

It hits your foot's trigger points perfectly! The biomechanical function of the **TMX[®] META** by **BlackBoard[®]** makes your foot suppler. The **TMX[®] META** is 60% mobilizer and 40% trigger point massager.

GENERAL USAGE

USAGE POSSIBILITIES

You can place the **TMX[®] META** vertically or horizontally under your foot.

VARIATIONS

By turning the **TMX[®] META**, you will be able to perfectly hit the point of your pain. **IMPORTANT:** The smallest of adjustments, or corrections, of the **TMX[®] META** can prompt fundamental improvements in the treatment. You can find many effective ways to use your **META** and advice on numerous syndromes in video tutorials on our website.

REACTIONS

Reactions to your personalized treatment are normal. Sore muscles, redness of the skin and localized pain are reactions in the tissue that may occur after usage. About 10-48 hours after using the **TMX[®] META** you might feel worse at first. Please be aware that this is only a temporary reaction! Make sure that you use your device exactly as shown in the images. Pay attention to your breathing at all times. Don't hold your breath. Bear in mind that the first 2-3 minutes could be painful. However, your muscle will adapt to the irritation, and then you'll feel relief. Once you've mastered an exercise, take it to the next level. If the intensity and duration of your pain increases, stop the exercise.

RELEASES MUSCLES IN 3 SIMPLE STEPS

1 TRIGGER 2 MOBILIZE 3 STRETCH

TRIGGER
Place your foot, standing, on the **TMX[®] META** to find the tensions or trigger points.

MOBILIZE
Stretch and flex your toes while trying to shift your body's center of gravity back and forth over your knee joint. This will create enough pressure and movement at the same time.

STRETCH
At the end, stretch out the soles of your feet and calf muscles.

DURATION OF USE
You can complete this **TMX[®] 3-step program** in about 5-10 minutes, as often as you want.

INDICATIONS

General foot issues, blocked feet, heel pain, foot tension, pain on the sole of the foot (*plantar fasciitis*), midfoot pain (*metatarsal pain*), heel spurs, hallux valgus, hallux rigidus (*Osteoarthritis* of the big toe in the metatarsophalangeal joint), pain on the inner part of the foot, pain on the outer part of the foot, bad foot statics, foot spasms (neurological), foot sole sensory impairment, pain on the ball of the foot, diabetic foot.

CONTRAINDICATIONS

Acute signs of inflammation such as swelling, redness, warmth, severe pain with inability to bear weight on the foot, acute trauma, acute supination trauma, ruptures of the ankle joint, fractures, edema, Charcot arthropathy, bone softening, osteoporosis, inflammatory impingement.

PAIN AROUND THE OUTER EDGE OF THE FOOT (E.G. AFTER SUPINATION TRAUMA)

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

For best results, move the **TMX[®] META** about 1-2 cm in every direction over the painful area. See if you can tilt your foot gently left and right (*Pronation/Supination*). Try placing the **TMX[®] META** wedge between your 4th and 5th metatarsal bones.

PAIN AROUND THE INNER EDGE OF THE FOOT (E.G. AFTER RUNNING / WRONG RUNNING SHOE)

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

For best results, move the **TMX[®] META** about 1-2 cm in every direction over the painful area. See if you can tilt your foot gently left and right (*Pronation/Supination*). Try placing the **TMX[®] META** wedge between your 1st and 2nd midfoot bones.

PAIN AROUND THE MIDFOOT (METATARSALGIA - SUPREME DISCIPLINE)

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

For best results, move the **TMX[®] META** about 1-2 cm in every direction over the painful area. See if you can tilt your foot gently left and right (*Pronation/Supination*). Move the **TMX[®] META** between each and every metatarsal bone, so you can loosen up the area optimally. It might hurt more between the individual metatarsal sections. Don't worry, that's how it should be. First it hurts, then it feels good!

HALLUX VALGUS (BIG TOEBENT OUTWARDS)

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

By moving the **TMX[®] META** at an angle, you can reach the oblique muscles that support the metatarsophalangeal joint. Here you can massage the trigger point in the most important muscle (*M. flexor hallucis brevis*). If this muscle is overstrained, the metatarsophalangeal joint will become unstable, which eventually leads to degeneration.

USAGE HALLUX RIGIDUS (METATARSOPHALANGEAL JOINT OSTEOARTHRITIS)

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

By moving the **TMX[®] META** at an angle, you can reach the oblique muscles that support your metatarsophalangeal joint. Here you can massage the trigger point in the most important muscle (*M. flexor hallucis brevis*). If this muscle is overstrained, the metatarsophalangeal joint will become unstable, which eventually leads to degeneration.

USAGE SCREWFOOT / FLAT- FEET / SPLAYFOOT / HIGH ARCHES

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

You can apply the **TMX[®] META** anywhere you feel pain - turn or move the **META** as desired. In congenital deformities, the foot muscles are considerably overworked. Give yourself enough time to alleviate points of tension little by little.

PAIN BALL OF FOOT („FOOT OPENER“)

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

You can apply the **TMX[®] META** anywhere you feel pain - turn or move the **META** as desired. Correct usage of the "foot opener" will result in unrivaled gait experiences. Give yourself enough time to work through the entire ball of your foot.

USAGE TAILOR'S BUNIONS (DIGITUS QUINTUS VARUS / PROMINENCE OF THE SMALL TOE)

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

You can apply the **TMX[®] META** anywhere you feel pain - turn or move the **META** as desired. Mobilize all of your metatarsal bones, dedicating more time to the area between the 4th and 5th bone. It might hurt more here. Don't worry, it's supposed to hurt more at first. First it hurts - then it feels good!

OSTEOARTHRITIS INDIVIDUAL TARSAL BONES

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

You can apply the **TMX[®] META** anywhere you feel pain - turn or move the **META** as desired. Start with the **TMX[®] META** around the area in front of your heel. This is where you can mobilize your tarsal bones in various ways.

PAIN IN THE HEEL / HEEL SPURS

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

You can apply the **TMX[®] META** anywhere you feel pain - turn or move the **META** as desired. Do not place the **TMX[®] META** only on your heel. Be sure to position it in front of the heel as well since the surrounding tissue also needs to be treated.

USAGE FEELING OF BEING BLOCKED

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

You can apply the **TMX[®] META** anywhere you feel pain - turn or move the **META** as desired. Don't be afraid if you hear a cracking sound. It is a completely normal occurrence when you mobilize the foot.

USAGE FOOT SOLE PAIN (PLANTAR FASCIITIS)

Usage according to the **TMX[®] 3-step program**

SUGGESTIONS

You can apply the **TMX[®] META** anywhere you feel pain - turn or move the **META** as desired. Start right in front of the heel, moving the **META** forward little by little towards the ball of the foot. This is an excellent way to work on your plantar tendon.

TMX[®] CONNECT

THE WAY TO THE TMX[®] WEBAPP

Sign up for our free **TMX[®] webapp** and receive a wealth of accompanying application tutorials and in-depth trigger knowledge: www.tmx-connect.de